

TECHNICAL DATA SHEET



Product Name and Code:

Type 1 Chemical Indicator for Hydrogen Peroxide – **LCH41**

Product Description:







LISTER Type 1 Chemical Indicator for VH_2O_2 LCH41 has an indicator ink bar printed on a white plastic strip. The chemical indicator bar turns from green to yellow after exposure to vaporized hydrogen peroxide.



Indications for use:

Use the LISTER Type 1 Chemical Indicator for VH_2O_2 LCH41 as an internal pack process indicator to verify exposure to vaporized hydrogen peroxide in the sterilizers. The chemical indicator bar turns from green to yellow after exposure to vaporized hydrogen peroxide.

Performance characteristics:

Reference Color Change Conditions	Illustration
2.3 mg/l VH_2O_2 , 3 mins	 UNPROCESSED
	 FAILURES One or more sterilization parameters were not achieved. The materials cannot be considered sterile and must be processed again. The main failure causes are: Inadequate cycle duration Improper cycle selection Overloading / incorrect loading of the sterilizer Sterilizer malfunction Blocking indicator with large metal mass
	 PROCESSED AND SATISFACTORY
	
	
	

Technical information:

Weight	117-121g / bag
Period of validity	36 months
Dimensions	1.7 cm×10 cm

Compliance:

ISO 11140-1:2014 Sterilization of health care products-Chemical indicators-Part 1: General requirements.

Manufactured by: LISTER BIOMEDICAL CO.,LTD.
Building 2, No.72 Shengli Road, Yixiu, Anqing, Anhui, China.
Tel. +86 556 5689070
E-mail: info@listerbiomed.com -www.listerbiomed.com

Package:

200PCS/PACK

Version: 20241209